

Product No. 18186

# TMB Solution for Western Blotting

TMB Solution for Western Blotting is a horseradish peroxidase (HRP) detection reagent for western blotting or immunoblotting.

## Features

- High detection sensitivity similar to general chemiluminescence substrate without any equipment.
- Produce visible dark blue bands (or dots).
- Ready-to-use
- No carcinogen such as 3, 3'-Diaminobenzidine or o-Dianisidine is used.

## Composition

TMB, Hydrogen Peroxide, Precipitate agent, Stabilizer

## Protocol

1. After the blotting with HRP-conjugated antibody, wash the membrane thoroughly with TBS-T.
2. Place the membrane in a plastic tray with the protein side up and add TMB solution to cover the membrane completely. (Required volume is 0.1 mL per cm<sup>2</sup> of membrane.)
3. Incubate for 5 - 15 minutes without shaking<sup>\*1</sup>.
4. When the expected band appears on the membrane, take the membrane out of the solution, and wash it with distilled water<sup>\*2</sup> several times to stop reaction.

<sup>\*1</sup> Shaking may decrease the sensitivity.

<sup>\*2</sup> Be sure to use distilled water for washing. Washing with other solution, e.g., TBS-T, decrease the sensitivity.

## Caution

- To store the membrane, put the dried membrane in a plastic bag and protect from light. Although the membrane is stored correctly, the band will fade gradually. It is recommended to scan the image within one week.
- In case of high background, the concentration of antibodies may be too high, or the washing step may be insufficient. Optimize the concentration of antibodies, or increase the number washing step.
- In case of low signal, optimize the concentration of antibodies or the amount of loading protein.
- This product is only for immunoblotting, and not for immunohistochemistry or ELISA.

## Storage

Refrigerator

## Expiration

Expiration date is stated on the product label.

## Packing

200 mL (Product No.18186-24)