

Human Pluripotent Stem Cell Passaging Protocol

This protocol is optimized for BD standard 6-well flat-bottomed tissue culture dishes. However, the protocol can be adapted to tissue culture vessels of different sizes and configurations.

Cell Passaging Steps

1. Pre-warm the required volume of cell culture medium containing 10 μ M Y-27632 (final concentration) to 37°C.
The typical volume is 2 ml for every well that the cells will be transferred to + 1 ml for each well that needs to be harvested + extra 1 - 2 ml.
Examples:
 1. When you pass the cells from 1 well to 6 wells:
-> $[(2 \text{ ml} \times 6 \text{ wells}) + 1 \text{ ml for harvesting} + \text{extra } 1 - 2 \text{ ml}] = 14 - 15 \text{ ml medium}$
 2. When you pass the cells from 2 wells to 8 wells:
-> $[(2 \text{ ml} \times 8 \text{ wells}) + 2 \text{ ml for harvesting} + \text{extra } 1 - 2 \text{ ml}] = 19 - 20 \text{ ml medium}$
2. Remove the spent culture medium.
3. Wash the adherent cells gently with 2 ml of D-PBS(-) (Thermo Fisher Scientific) in each well.
4. Add 600 μ l 0.5x TrypLE Select (Thermo Fisher Scientific) diluted in 0.5 mM EDTA (Sigma-Aldrich)/D-PBS(-) to each well. **See below on how to prepare 0.5 mM EDTA/D-PBS(-).*
5. Incubate the plate(s) in an incubator for 10 minutes at 37°C, 5% CO₂.
6. Aspirate the supernatant (TrypLE Select + EDTA) from each well. Wash the adherent cells remaining on the plate gently with 2 ml D-PBS(-).
7. Add 1 ml of the pre-warmed medium (37°C) with the 10 μ M Y-27632 to each well.
8. Harvest the cells using a cell scraper, followed by gentle pipetting to generate a single cell suspension. **Being gentle with the scraper is critical for retaining high cell viability.**
9. Stain the cells with 10 μ l of Trypan blue staining solution (Thermo Fisher Scientific). Determine the number of viable cells either by counting manually with a haemocytometer or by using an automated cell counter.
10. Adjust the concentration of the cell suspension with the pre-warmed medium containing 10 μ M Y-27632. The recommended cell density is 2.0 - 5.0 $\times 10^4$ cells/well in a 6-well plate, depending on the hPSC line used (Miyazaki et al., 2012;

Miyazaki et al., 2017). The number of cells will increase 100-fold after 6 - 7 days in culture.

11. Add an adequate volume of iMatrix-511. The recommended iMatrix-511 concentration is 0.25 $\mu\text{g}/\text{cm}^2$ (4.8 μl of 500 $\mu\text{g}/\text{ml}$ iMatrix-511/well).
12. Transfer 2 ml of the cell suspension with the iMatrix-511 to each well.
13. Gently rock the plate back-and-forth and side-to-side to disperse the cells across the surface.
14. Incubate the cells for 24 hrs before exchanging the plating medium for fresh medium (without Y-27632). Replace medium on day 3 and day 5, then everyday until the next passage (usually at day 6 or day 7).

Recommended Volume of Reagents

	6 well	12 well	24 well	48 well
Approximate Area (cm^2)	9.6	4.0	2.0	1.0
0.25 $\mu\text{g}/\text{cm}^2$ iMatrix-511 (μl)	5	2	1	0.5
Medium volume (ml)	2	1	0.5	0.25
TrypLE Select/EDTA/D-PBS(-) (μl)	500	250	125	65

*Preparation of 0.5 mM EDTA/D-PBS(-)

1. First, prepare the 0.5 M EDTA/D-PBS(-) stock solution by adding 18.6 g EDTA-2Na (Sigma-Aldrich) to 100 ml D-PBS(-). Mix well until EDTA-2Na is completely dissolved. Adjust the pH to 8.0. Store at room temperature, in the dark, until needed.
2. To prepare the 0.5 mM EDTA/D-PBS(-) working solution, use the concentrated solution in step #1 and make a 1,000-fold dilution by adding D-PBS(-). Store it at 4°C until use. This solution is good for two weeks, any remaining solution after two weeks should be safely disposed.



Nacalai USA, Inc.

6625 Top Gun Street, Suite 107,
San Diego, CA 92121
TEL: 858 404 0403
Email: info@nacalaiusa.com
Website: www.nacalaiusa.com