


HistoVT One

INTRODUCTION

HistoVT One is an antigen retrieval solution for immunohistochemistry and in situ hybridization. This product can unmask antigenic sites without damage to antigen from formalin-fixed, frozen or paraffin-embedded tissue sections.


- Enhancing antigen-antibody reaction
- Applicable to frozen or paraffin-embedded tissue sections
- High reproducibility

APPLICATION DATA 1 Formalin-fixed, paraffin-embedded tissue sections



Blue: PCNA (Immunohistochemistry), Green: Nuclear staining (Cytogreen staining), Red: Telomere DNA

APPLICATION DATA 2 Frozen section



Data courtesy of RIKEN Brain Science Institute, Brain Development Research Group.

ORDERING INFORMATION

Item Name	Storage	Cat. No.	Qty	Price
HistoVT One (10x, pH 7.0)	RT	06380-05	500 ml	\$200