



Microplate Clear Seal for Real-Time PCR

\$1.35 Unit price
per sheet

Unique technology Makes it easy to peel and stick



With low autofluorescence and it can
be used for PCR and real-time PCR

The product is recommended for those who:

- ✓ Want to reduce the cost of disposables
- ✓ Want to use high quality products
- ✓ Want to use easy-to-apply products

Manufacturing in super clean room ISO class1



Cat. No.	PKG Size	Price (US\$)
R-CSC004-01	100 sheets	135.00
R-CSC004-02	5 x 100 sheets	608.00

User voice (Stem Cell Biologist)

In real-time PCR experiment, there is a concern that performance (cost-effectiveness) decreases if the cost of consumable materials is reduced. However, this clear seal is cheaper than the conventional one and the data accuracy is satisfactory.

This type of experiment is always large in quantity, so the total cost has dropped significantly. In addition, since it is stably supplied through strict quality control in production, it does not waste valuable experiments and samples under any circumstances.

Based on this background, this seal has now become the standard for this experimental protocol in our laboratory and has become an indispensable consumable article.

Applications: for PCR, real-time PCR
(Pressure bonding type)

Materials: PET

Size: 138 mm x 78 mm
(without tab: 119 mm x 78 mm)

Temp: -20 to 110 degree C