

Maximum flexibility for designing your own microplate

CS CRIE's

Bottomless Microplate

Attach any materials to the bottom to make original microplate!

Bottomless microplate is a first choice to make an original microplate with different material of bottom.

Our quality bottomless microplate with maximum flexibility accelerates the design of original microplate that meets your experimental design.

Available in 4 Formats with different tapes




- 6, 24, 96 and 384-well microplate
- No tape, Acrylic adhesive tape or Silicone adhesive tape on the bottom
- Material : Polystyrene



White: 96 & 384-well



with adhesive tape



1536-well

Coming Soon!

Product and Price List

CSCRIE Bottomless Microplate, No Tape (5 plates)		CSCRIE Bottomless Microplate, with Acrylic adhesive Tape (5 plates)		CSCRIE Bottomless Microplate, with Silicone adhesive Tape (5 plates)	
Product Code	Well Type	Product Code	Well Type	Product Code	Well Type
CSCSM0061	6-well	CSCSM0062	6-well	CSCSM0063	6-well
CSCSM0241	24-well	CSCSM0242	24-well	CSCSM0243	24-well
CSCSM0961	96-well	CSCSM0962	96-well	CSCSM0963	96-well
CSCSM3841	384-well	CSCSM3842	384-well	CSCSM3843	384-well

50 plates and more are available, please contact us.



Custom Services!

We have extensive expertise and technologies for bonding and creating customized products!

Bonding own materials (e.g. Glass, PET / PDMS / COP film and etc.) firmly to the bottomless plate needs enough experience and technique. Insufficient bonding is the main factor of inter-well leaking and contamination.

CS TEC, the mother company of CS CRIE, has over 23 years of experience for manufacturing of a variety of microplates. Moreover, our all latest equipment for punching, laser micro-processing and bonding enable us to manufacture quality products and customized products that meet customer's various requirements

Not only microplates, but specially processed microfluidic chips are also offered.


Contact us now!



Examples of Custom Products



Micro-processed bottom in the culture plate



Microplate with fluidic channels



Zero-flat bottom and perfusion channel



Original cell culture insert



Microfluidic plate



3-part plate



Learn more : www.cscrie-eng.com/contents/bottomless_microplate



Manufacturer : CS CRIE CO., LTD.

57-1 Ogahara, Makishimacho, Uji-shi, Kyoto 611-0041, Japan

TEL.+81-774-29-9037

E-mail: yagi@cstec-jp.com

Website : <https://www.cscrie-eng.com>

Feel free to contact us about any question, quote and sample request and customization project!