


Human polyclonal antibody prepared from serum of SARS-CoV-2 convalescent




Here, EVEC Inc. prepared immunoglobulin fraction from serum of SARS-CoV-2 convalescent.

The product is categorized into human polyclonal antibody, and DOES NOT include any infectious materials such as virus particles because of filtration by 0.22 μm membrane.


Of course, as the product is consisted of total immunoglobulin that purified by Protein A affinity column, various Igs are contained, including anti-Spike protein and anti-Nucleocapsid etc..

Please use this product for research and elucidation the pathology of COVID-19 symptom.



Western blot

Western blot analyses of this product against SARS-CoV-2 S-protein (S) or N-protein (N)



ELISA

Standard curve of this product against SARS-CoV-2 Spike protein

Product Name	Type	Application	Catalog Number	PKG Size
Total Ig of SARS-CoV-2 convalescent	Human-Poly	WB, ELISA	EVCI-001-0050	50 μg
			EVCI-001-0100	100 μg
			EVCI-001-0500	500 μg

EVEC, Inc. Kobe Laboratory

Kitakan 404 Kobe International Business Center
5-5-2, Minatojima minamimachi, Chuo-ku,
Kobe-shi, Kobe Hyogo, 650-0047 JAPAN
E-mail : kobe@evec.jp Web : <https://www.evec.jp/>

Distributor

NACALAI USA, INC.

10225 Barnes Canyon Rd., Suite A103
San Diego, CA 92121
E-mail : info@nacalaiusa.com
Web : www.nacalaiusa.com