


EZSPHERE Variety Pack Microplate 6-well EZ-TRY

TCI-4810-EZ-TRY-SP-N

5 plates per box

EZSPHERE Plate: 6 types of micro-well are created on the surface of 6-well plate.

Micro-well specification: Please see the details below.



Enable efficient formation of uniformly-sized spheroids in a large quantity. Low adhesion is performed on the surface of micro-fabricated wells. EZSPHERE series have 6 types of micro-fabricated wells.

Type	Dia. X Depth (μm)	Micro-well numbers
900	400-600 x 100-300	2,700
901	250 x 100	10,600
902	450 x 200	2,700
903	950 x 400	1,300
904	1,000 x 350	700
905	1,600 x 800	260

* Micro-well size and numbers are approximate value.


EZSPHERE Variety Pack Microplate 24-well EZ-TRY

TCI-4820-EZ-TRY-SP

5 plates per box

EZSPHERE Plate: 6 types of micro-well are created on the surface of 24 well plate.

Micro-well specification: Please see the details below.



Enable efficient formation of uniformly-sized spheroids in a large quantity. Low adhesion is performed on the surface of micro-fabricated wells. EZSPHERE series have 6 types of micro-fabricated wells.

Type	Dia. X Depth (μm)	Micro-well numbers
900	400-600 x 100-300	470
901	250 x 100	1,840
902	450 x 200	470
903	950 x 400	210
904	1,000 x 350	120
905	1,600 x 800	45

* Micro-well size and numbers are approximate value.