

Temperature Controlled Package TACPack-20 Series

Thermo
Accumulation
Carriage
Pack

PRODUCED BY TAMAIKASEI Co., Ltd
<http://www.tamai-kasei.co.jp>
made in JAPAN

TACPack-20 Series are parcel sized temperature controlled packages that maintain the load temperature below -20°C. Simple to pack and with a selection of sizes, it provides easy solution to keep your products within the right temperature range.


Product Specification

Product	Package Size (mm)	Payload Size (mm)	Payload Capacity (ℓ)	Weight (kg)	Material	Outer	PCM	Q'ty
TACPack-20SS	330x220x250	190x20x135	0.5	1.6	EPS	Carton	Cold Ice CV 30 500g	2
TACPack-20S4	520x320x395	230x180x80	3.3	8.2	EPS	Carton	Cold Ice CV 30 1300g	4
TACPack-20S5	520x320x395	230x180x45	1.8	9.7	EPS	Carton	Cold Ice CV 30 1300g	5
TACPack-20H8	480x455x535	180x220x195	7.7	15.7	EPS	Carton	Cold Ice CV 30 1300g	8
TACPack-20H9	480x455x535	180x220x155	6.1	17.2	EPS	Carton	Cold Ice CV 30 1300g	9
TACPack-20L	640x450x500	355x215x163	12.4	25.7	EPS	Carton	Cold Ice CV 30 1300g	14

Duration

Product	Temperature Range	35°C	25°C	-10°C
TACPack-20SS	below -20°C	13hr	16hr	62hr
TACPack-20S4		20hr	37hr	163hr
TACPack-20S5		49hr	59hr	168hr
TACPack-20H8		33hr	41hr	168hr
TACPack-20H9		52hr	62hr	168hr
TACPack-20L		59hr	66hr	168hr

※The data above does not guarantee its performance.



TAMAIKASEI Co., Ltd.

Web: <http://www.tamai-kasei.co.jp> Email: mailto@tamai-kasei.co.jp
3-524-9 Zenibako Otaru, Hokkaido, JAPAN 〒047-0261
TEL : +81-134-62-1100 FAX : +81-134-62-0777